




FAST RESULTS

***FLEXIBLE
THROUGHPUT***

***FULL
CUSTOMIZATION***

- ⚙ Fully automated from sample to result
- ⚙ Versatile and fast diagnostics platform for 1 – 16 samples
- ⚙ High multiplex qPCR capabilities
- ⚙ Full bi-directional LIS connectivity options



WHAT MAKES IT UNIQUE & EXCEPTIONAL!

Patented Spin-Tip Technology, enabling fast protocols (< 30 minutes) with optimal results

In-Tip magnetic bead extraction technology for enhanced efficiency

High Multiplex 5+ channel qPCR unit for up to 32 reactions (split-sample option)

Integrated PCR set-up for any lyophilized qPCR solution

Single-strip extraction reagents, optimized for selected applications

Primary Sample Handling Module (incl. LLD and barcode reading) optional



 **ZINEXTS** NeXgen 16HP
Zinexts Life Science Corp.

COMING SOON !

For further inquiries, please contact sales@zinexts.com

Contact us

www.zinexts.com
sales@zinexts.com
+886 2246-3579



16F., No. 93, Sec. 1, Xintai 5th Rd., Xizhi Dist., New Taipei City, Taiwan

ZFL-NeXgen-0123-V1