

UniColl

Sample Collection Tube

ONE Fits All



The savior to your troublesome and complex pretreatment process

Sample treatment is one of the most important parts of the genomic process, and effective sample preparation is essential for achieving good analytical results.

We provide the ideal sample preparation method that is fast, accurate, precise and keeps sample integrity.



Key Features

- ④ Caps distinguished by colors contain optimized buffers, designed for your different sample pretreatment.
- ④ Universal size for centrifuge & vortex.
- ④ Simple, fast, and contact free. Provides minimal hands on time for sample preparation and collection.
- ④ Great sample quality for extraction.

How It Works



Sample Collection


Convenient and fast collection of samples with different collecting tools, depending on sample types.



Sample Pretreatment

Simply twist the cap to mix buffer with sample.

With tube size designed for universal centrifuge & vortex, no sample transferring is required for pretreatment.



Sample Preparation

Special buffers custom made for different sample types are stored in the twistable cap.

Different colored caps provide quick and easy identification of tube purpose and buffers.



Ready for Extraction

Contact-free sample load with just a few drops by pulling off the bottom cap and removing the lock.

What we have

Stool



TB



Tissue, Cell



Urine



Saliva



With more to come

Contact us

www.zinexts.com

Sales@zinexts.com

+886 2246-3579

16F., No. 93, Sec. 1, Xintai 5th Rd., Xizhi Dist., New Taipei City, Taiwan

ZFL- SampleTube -0919-V3