


AUTOMATED TOTAL SOLUTION

for Molecular Diagnostics



**FAST
RESULTS**

**FLEXIBLE
THROUGHPUT**

**FULL
CUSTOMIZATION**

- Fully automated from sample to result
- Versatile and fast diagnostic platform for 1 – 16 samples
- High multiplex qPCR capabilities
- For lyophilized and liquid qPCR assays



WHAT MAKES IT UNIQUE & EXCEPTIONAL!

- Patented Spin-Tip Technology, enabling fast protocols (< 30 minutes) with optimal results
- In-Tip magnetic bead extraction technology for enhanced efficiency
- High Multiplex 5+ channel qPCR unit for up to 32 reactions (split-sample option)
- Integrated PCR set-up for any qPCR solution (lyophilized & liquid form)
- Primary Sample Handling Module (incl. LLD) optional
- Integrated 2D Barcode Reader for easy ID scanning
- Multiple connectivity options, incl. bi-directional LIMS

